



World-Best Ultrasound Paper,
 Fitting Perfectly to World-Renowned Thermal Printers, Sony & Mitsubishi

SUPER ULSTAR Ultrasound Paper



FDA Registered
 Certified to ISO 9001 / ISO 14001 / ISO 13485



DURICO DURICO C&T, Inc.
<http://www.durico.co.kr>



Why pay more for the same quality?

Compatible Papers for Mitsubishi Printer

- ▶ No need to make any adjustment on the printer for optimum printed images



Paper Grade, Size and Compatibility

Grade	Size	DURICO	SONY	MITSUBISHI
High Glossy	110mmX18m	ULSTAR-1100HG	UPP-110HG	K91HG/KP91HG
	84mmX12.5m	ULSTAR-840HG NEW	UPP-84HG	
High Density	110mmX20m	ULSTAR-1100HD	UPP-110HD	
	110mmX20m	ULSTAR-1100HD mibi		K65HM/KP65HM
	210mmX25m	ULSTAR-2100HD	UPP-210HD	
High Density Matt	110mmX20m	ULSTAR-1100HD Matt		
Standard	110mmX20m	ULSTAR-1100S	UPP-110S	
	110mmX20m	ULSTAR-1100S mibi		K61S/KP61S/KP61B
	84mmX13.5m	ULSTAR-840S NEW	UPP-84S	

Packing Information

- ▶ Economic Purchase for smaller Packing Unit (5 rolls in an inner box)



- ULSTAR-1100 series:** Each roll in a pouch
 5 rolls in an inner box
 10 inner boxes in a master carton (50 rolls)
- ULSTAR-2100HD:** Each roll in a pouch
 5 rolls in an inner box
 5 inner boxes in a master carton (25 rolls)
- ULSTAR- 840 series:** Each roll in a pouch
 10 rolls in an inner box
 10 inner boxes in a master carton (100 rolls)

DURICO Durico Coating & Technology, Inc.

Seoul Office Banpo-daero 26-gil 21, Seocho-gu(Seocho-dong, Sausohun Bldg. 2F), Seoul 06648, Korea Phone: 82-2-525-8405 Fax: 82-2-525-7461 E-mail: info@durico.co.kr	H/Q and Plant Oedap 6-gil 33, Sangju-si, Gyeongbuk 37240, Korea E-mail: info@durico.co.kr	U.S. Office Durico Imaging, Inc. 18553 S.Dominguez Hills Dr. Rancho Dominguez, CA 90220, U.S.A Office: 714-521-7000 Direct: 714-868-1945 E-mail: sales@duricoimaging.com	Europe Office Durico Imaging s.a.r.l. 158 Rue Diderot, 93500 Pantin, FRANCE Office: +33(0)9 82 56 98 59, Fax: +33(0)9 82 59 38 59, E-mail: euro@duricoimaging.com
--	--	---	--


Direct Thermal Film for Video Printer

Durico is specialized in Self-Manufacturing Top Quality Thermal Print Media for printing Ultrasound Scanned Images.



► Totally different from other Inferior Compatible Papers

► Highly Protective UV-Cured Top Layer




Long Lasting Image Durability for High Storage Stability

► High Humidity and Heat Resistance for Long Storage


Humidity resistance 40°C, 90%, 7Days		Heat resistance 50°C, Dry, 7Days	
ULSTAR	Other Compatible	ULSTAR	Other Compatible

► Strong Resistance to Water, Finger Prints, and Ultrasound Gel Drops



Precise, High Resolution Image Reproduction

► First-class Colour Intensity from Highest Black Colour Density




► Outstanding, Bright, Sharp & Detailed Medical Image Printing available from Precise Gradation of Grey-scale

► High Glossy, Brilliant Photorealistic Reproduction



Trouble-free Printing & Long-Term Use of Printer

► Easy Tearing along the Cross Direction



► Smooth Head Matching Performance for Continuous, Peace of Mind Printing



► Minimal Curling for Flat Photo Media & Easy Filing



► Built-in Anti-electrostatic Layer Preventing Thermal Head Damage

Durico

Other Compatible