



Why pay more for the same quality?

World-Best Ultrasound Paper

Perfectly Compatible with Sony & Mitsubishi Video Printers

SUPER ULSTAR

ULTRASOUND PAPER



DURICO DURICO C&T, Inc.

www.durico.co.kr


Direct Thermal Film for Video Printer

Durico is specialized in Self-Manufacturing Top Quality Thermal Print Media for printing Ultrasound Scanned Images.



▶ Totally different from other Inferior Compatible Papers

▶ Highly Protective UV-Cured Top Layer



Long Lasting Image Durability for High Storage Stability

▶ High Humidity and Heat Resistance for Long Storage


Humidity resistance 40°C, 90%, 7Days		Heat resistance 50°C, Dry, 7Days	
ULSTAR	Other Compatible	ULSTAR	Other Compatible

▶ Strong Resistance to Water, Finger Prints, and Ultrasound Gel Drops




Precise, High Resolution Image Reproduction

- ▶ First-class Colour Intensity from Highest Black Colour Density




- ▶ Outstanding, Bright, Sharp & Detailed Medical Image Printing available from Precise Gradation of Grey-scale

- ▶ High Glossy, Brilliant Photorealistic Reproduction



Trouble-free Printing & Long-Term Use of Printer

- ▶ Easy Tearing along the Cross Direction



- ▶ Smooth Head Matching Performance for Continuous, Peace of Mind Printing



- ▶ Minimal Curling for Flat Photo Media & Easy Filing



- ▶ Built-in Anti-electrostatic Layer Preventing Thermal Head Damage



Products

■ ULSTAR-1100 Series for SONY



- ▶ ULSTAR-1100HG | 110mm × 18m
- ▶ ULSTAR-1100HD | 110mm × 20m
- ▶ ULSTAR-1100S | 110mm × 20m
- ▶ ULSTAR-1100HD Matt | 110mm × 20m

Compatibility - UP-895 / 897 / 898 Series

Packaging - 5 rolls/box 50 rolls/carton

■ ULSTAR-1100 Series for MITSUBISHI



- ▶ ULSTAR-1100HG | 110mm × 18m
- ▶ ULSTAR-1100HD mibi | 110mm × 20m
- ▶ ULSTAR-1100S mibi | 110mm × 20m

Compatibility - P93 / P95

Packaging - 5 rolls/box 50 rolls/carton

■ ULSTAR-840 Series for SONY



- ▶ ULSTAR-840HG | 84mm × 12.5m
- ▶ ULSTAR-840S | 84mm × 13.5m

Compatibility - UP-D711MD

Packaging - 10 rolls/box 100 rolls/carton


■ ULSTAR-2100 Series for SONY



- ▶ ULSTAR-2100HD | 210mm × 25m

Compatibility - UP-970 / 971 / 990 / 991AD

Packaging - 5 rolls/box 25 rolls/carton



■ Seoul Office

Banpo-daero 26-gil 21,
Seocho-gu(Seocho-dong, Sausohun Bldg. 2F),
Seoul 06648, Korea
Phone: 82-2-525-8405 Fax: 82-2-525-7461
E-mail: info@durico.co.kr

■ H/Q and Plant

Oedap 6-gil 33, Sangju-si,
Gyeongbuk 37240, Korea
E-mail: ulstar.sales@durico.co.kr

■ U.S. Office

Durico Imaging, Inc
6905 Oslo Circle, Suite G
Buena Park, CA 90621, U.S.A.
Office : 714-521-7000 Direct: 714-868-1945
E-mail: sales@duricoimaging.com

■ Europe Office

Durico Imaging BVBA
Villastraat 2C, 1830 Machelen, Belgium
Office: +32 470 632 965,
Fax: +32 (0)2 791 9514,
E-mail: euro@duricoimaging.com